



User Manual A9

Please Read the Manual before Use

	Table of contents	1
	Product description	2
	Parts name	3
	Introduction for wired/wireless mode	4
	How to charge the battery	4
	How to apply the needle cartridge	5
	Adjustment of needle length	6
	Speed display and adjustment	6
	Needle cartridge introduction	7
	Reference chart	8
	Specification	8
	Precaution	9
	Warranty	10

Product description 2

The device consists of an electric microneedling pen and disposable needle cartridges. This manual will explain how to use each component in detail.

The device is used to improve the appearance of the skin and reduce signs of aging. It is intended to stimulate the skin with micro-needles for collagen and elastin production. The micro-channels created on surficial skin by device also help with the absorption of skin serums.

There are 7 different types of needle cartridges designed for the device. These cartridges can be used for micro-needling treatment and nano-needling treatment.

Parts name 3

Introduction for wired/wireless mode 4

Wired Model: Dr.pen A9-C
This model is able to work only if it is connected to the AC outlet. Connect the device to the adapter with the power cable, and plug the adapter into a live AC outlet to work.

Wireless Model: Dr.pen A9-W
This model is powered by a built-in battery. It can work wirelessly. The power cable and the adapter are used to charge the device when it is out of battery.

How to charge the battery

Connect the device to the adapter with the power cable, and plug the adapter into a live AC outlet to charge the device. Charging will be complete within 1 hour.

Charging Status



How to apply the needle cartridge 5

Preparation

1. Take out needle cartridge from the sterilized packaging. (Fig1)

Adjustment of needle length 6

Use suitable needle size according to different treating area and skin problem.

Adjust the depth of needle by turning the scale adjustment ring. The depth of needle should be adjusted after starting the pen. The scale number aligning with the scale point shows the current needle depth. (Fig3)

Needle cartridge introduction 7

Available needle cartridges

Reference chart 8

REGION	NEEDLE DEPTH (MM)
Forehead	0.25 - 0.5
Between the Eyebrow	0.25 - 0.5
Nose	0.25
The Area around Eye	0.25
Check and Chin	0.5 - 1.0
Surgical and Scar Maximum	2.5

Precaution 9

Please be sure to read the instructions before use and follow as described.

Do not apply it to the person as following:

1. Who has acute disease or bloody disease over 38°C
2. Who has allergy and hemophilia
3. Who has melanosis of the face
4. Who has high internal heat
5. Who regularly visits hospital
6. Who takes medicines
7. Who is under 12 years old

Do not use it on the power other than AC 220V. It may cause fires or malfunction.

Do not use it near the place with high temperature, dust or humidity (Near the humidifier). It may cause malfunction.

Do not repair, disassemble or alter the controller and the body. It may cause malfunction.

Do not touch the device or detach the power cords with wet hands. It may cause electric shock.

Do not use the stripped cord or propped cord. It may cause fires(Be sure to replace it with another).

Do not use the propped cord which has crack or scars. It may cause scars on the skin.

When you don't use it for a long time, unplug it from the outlet.

Be sure to have an inspection for test. If any troubles are found, get it repaired.

No one other than the operator can touch it. It may cause accidents or malfunction.

~	Alternating current	SN	Serial No
---	Direct current	EC REP	EU Representative
□	Class I equipment	CE	CE mark
⚠	Attention, consult accompanying documents	Manufacturer	Manufacturer
📅	Date of manufacture	i	Consult instruction for use

Warranty 10

The warranty period of the device is 12 months. We offer repair service from the purchase date. If any of our products has damages or natural trouble, we provide replacements or repair service. The use of third-party cartridges can cause damage to your device. This could invalidate your warranty.

General maintenance

Make sure to inspect the device and adapter before use.

Keep the device with power off. Wipe its surface out with a disinfecting wipes.

Pay attention to clean areas such as crevices, seams, and areas around where the needle cartridge attaches to the device.

Keep it away from direct sunlight or high temperature and humidity.

Do not leave the measurement, repair or alter of the device to non-professional staff.

When you use the device which has not been used for a long time, make sure to check if it works properly with safety.

2. Find the "U" shape notch on the cartridge, align the notch with the dot on the bottom of the device, and slightly push the cartridge into the device. (Fig2)
3. Lock the needle cartridge to the device by turning it clockwise.
4. Turn the needle cartridge counterclockwise and pull it out when you would like to remove the cartridge.
5. Start the device only if the needle cartridge is assembled correctly. Press and hold the power button for 2-3 seconds to turn on/off the device.

Before assembling or disassembling the needle cartridge, make sure the U-shaped point of the needle cartridge is aligned with the point of the handpiece.

The needle cartridge is designed for single use.

Place the used needle cartridges in a trash container. If the needle cartridge is inadvertently contaminated during installation (i.e. dropped on the floor, package opened/damaged. The needle cartridge may be contaminated) Please discard it and get a new needle cartridge.

Speed display and adjustment 6

LED light speed display

The blue light turns on when the machine starts.

6 Speeds Adjustment

Press the button to change the speed. Turning on one more light means the speed increases 1 level.

Attention

The needle cartridge is designed for Dr.pen A9 only. The needle cartridge is disposable. Please sterilize the needle cartridge with 75% alcohol for few minutes before use. Please dispose the used needle cartridge to the trash container.

Specification 8

Product name: Auto-stamp Motorized Meso Machine
Model name: Dr.pen A9
Adapter
Input: 100-240V 50-60HZ
Output: 5V-1000mA
Package Weight: 420 g
Package Size: 169*123*56 mm

Packing List

Item	Quantity
Hand Piece	1 pc
Adapter	1 pc
Manual	1 pc
Power Cable	1 pc
Needle Cartridge	2 pcs

Precaution 9

Please be sure to read the instructions before use and follow as described.

Do not apply it to the person as following:

1. Who has acute disease or bloody disease over 38°C
2. Who has allergy and hemophilia
3. Who has melanosis of the face
4. Who has high internal heat
5. Who regularly visits hospital
6. Who takes medicines
7. Who is under 12 years old

Do not use it on the power other than AC 220V. It may cause fires or malfunction.

Do not use it near the place with high temperature, dust or humidity (Near the humidifier). It may cause malfunction.

Do not repair, disassemble or alter the controller and the body. It may cause malfunction.

Do not touch the device or detach the power cords with wet hands. It may cause electric shock.

Do not use the stripped cord or propped cord. It may cause fires(Be sure to replace it with another).

Do not use the propped cord which has crack or scars. It may cause scars on the skin.

When you don't use it for a long time, unplug it from the outlet.

Be sure to have an inspection for test. If any troubles are found, get it repaired.

No one other than the operator can touch it. It may cause accidents or malfunction.

~	Alternating current	SN	Serial No
---	Direct current	EC REP	EU Representative
□	Class I equipment	CE	CE mark
⚠	Attention, consult accompanying documents	Manufacturer	Manufacturer
📅	Date of manufacture	i	Consult instruction for use

Warranty 10

The warranty period of the device is 12 months. We offer repair service from the purchase date. If any of our products has damages or natural trouble, we provide replacements or repair service. The use of third-party cartridges can cause damage to your device. This could invalidate your warranty.

General maintenance

Make sure to inspect the device and adapter before use.

Keep the device with power off. Wipe its surface out with a disinfecting wipes.

Pay attention to clean areas such as crevices, seams, and areas around where the needle cartridge attaches to the device.

Keep it away from direct sunlight or high temperature and humidity.

Do not leave the measurement, repair or alter of the device to non-professional staff.

When you use the device which has not been used for a long time, make sure to check if it works properly with safety.