



ULTIMA A7

User Manual



Parts Name

How to Operate

- 1. Setting**
 -Adapter Mode
 (1) Plug in DC adaptor to the socket.
 (2) Plug in DC adaptor to AC outlet.
- 2. Preparation**
 (1) Take out needle cartridge in the sterilized packing.
 (2) Fitting round slots 1 into 2 of machine, push needle cartridge horizontally up to the end.
-
- (3) Lock the cartridge by turning it to the counter clockwise.
 (4) Adjustment of needle length should be done after turn on the machine.
 (5) If you want to remove the needle cartridge, turn it to the clockwise and push out horizontally, fitting round slots(slot1.slot2) together.

How to Operate

- 3. Power Switch**
-
- 4. Adjustment of speed**
-
- Rotate the speed ring to adjust the speed.
- The speed will be faster if you turn the ring clockwise.
 - The speed will be slower if you turn the ring counterclockwise.

How to Operate

- 5. Adjustment of needle length**
- You can adjust the length of needle by turning the adjustment ring clockwise or counter clockwise. Turning it clockwise, length of it will be longer. Turning it counter clockwise, length of it will be shorter.
-
- Reference Chart**
- | REGION | THIN SKIN | THICK SKIN |
|-----------------------|------------|-------------|
| Forehead | 0.25-0.5mm | 0.5-0.75mm |
| Between Eyebrow | 0.25-0.5mm | 0.5-1.0mm |
| Nose | 0.25mm | 0.5mm |
| Eye Area | 0.25mm | 0.25-0.5mm |
| Cheek Bone | 0.5mm | 0.5-1.0mm |
| Cheek and Chin | 0.5-1.0mm | 1.0-2.0mm |
| Upper Lip | 0.25mm | 0.25-0.75mm |
| Facial Scarring | 1.0-1.25mm | 1.5-2.0mm |
| Stretch Mark | 1.5-2.0mm | 1.5-2.0mm |
| Surgical & Burn Scars | 2.0-2.5mm | |

Specification

- 1. Object**
 Apparatus for applying a Micro-needle Therapy by vibration.
- 2. Performance**
 A device for puncturing skin-penetrating drugs or cosmetics supports with a micro-needle guided between extended and retracted positions.
- 3. Classifications**
 -Type of protection against electric shock: Class II Equipment
 -Degree of protection against electric shock: Floating Applied Part
 -Degree of protection against the ingress of water: Ordinary equipment(IPX0)
 -Mode of operation: Continuous operation
- 4. Specification**
 1)Product name: Auto-stamp Motorized Meso Machine
 2)Model name: Ultima A7
 3)Adapter: 5V, 1000mA
 4)Weight: Main body 70g
 5)Size: 160*23 (mm)

Precaution

	Please be sure to read instructions before use and follow as described.
	Do not apply it to the man as follow: 1. Who has acute disease or bloody disease over 38℃ 2. Who has allergy and hemophilia 3. Who has redness of the face 4. Who has high internal heat 5. Who regularly visits hospital 6. Who takes medicines 7. Who is under 12 years old
	Do not use it on the power other than AC 220V. It may cause fires or malfunction.
	Do not use it near the place with high temperature, dust or humidity (Near the humidifier).It may cause malfunction.
	Do not repair, disassemble or alter the controller and the body. It may cause malfunction.
	Do not touch the device or detach the power cords with wet hands. It may cause electric shock.
	Do not use the stripped cord or frayed cord. It may cause fires or malfunction.
	Do not use the frayed cord which has crack or scars. It may cause scars on the skin.
	When you don't use it for a long time, unplug it from the outlet.
	Be sure to have an inspection for test. If any troubles are found, get it repaired.
	No one other than the operator can touch it. It may cause accidents or malfunction.

	Alternating current		Serial No
	Direct current		EU Representative
	Class II equipment		CE mark
	Attention, consult accompanying documents		Manufacturer
	Date of manufacture		Floating applied part

Warranty

- We guarantee the products as follow.
- For troubles of the products in use, contact _____.
- We provide the after-sale service within 7days from the receipt of an application form.
- The period when the parts of derma pen is 6 months after stopping of the production of the parts.
- We offer free repair to our products within 600 hours basis(warranty) from the day they are purchased. If any of our products has damages or natural trouble, we offer repairs or replacements.
- If they are included in the following problems, you have to pay for the repair even during warranty period.
 - Troubles due to a user's careless treatment, voluntary repair or alteration
 - Trouble due to the repairs by a man other than the technician from a dealer shop or the customer center appointed by manufacturer
 - Troubles or damages due to shocks while in movement
 - Trouble due to the accessories or options not appointed by manufacturer
 - Trouble due to other natural disaster(fires, salt or water, etc).

Maintenance

- 1. General maintenance**
 -Keep it with its power off. Wipe its surface out with a dry cloth.
 -Keep it away from direct sunlight or high temperature and humidity.
 -Do not leave the measurement, repair or alter of the device to non-professional staff.
- 2. Routine inspection**
 -Make sure to inspect the device and adapter before use.
 -When you use the device which has not been used for a long time, make sure to check if it works properly with safety.
- Packing list**
- | | |
|------------|-------|
| Host | 1 pcs |
| Cartridges | 2 pcs |
| Manual | 1 pcs |
| Adapter | 1 pcs |

ULTIMA A7

User Manual